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****English - Dental Water Flosser Safety Warning****

****WARNING:** ** Risk of injury and infection. Please read the following instructions carefully before using this dental water flosser.

* ****Intended Use:** ** This product is intended for home oral hygiene use. Do not use for any other purpose. * ****Important Safety Instructions:** ** * Before each use, check that the device and its accessories are in good condition. Do not use if damaged. * Do not use this product if you have serious gum problems, excessive bleeding, or pain. Consult a dentist. * Never direct the water jet directly into the eyes. * Keep out of reach of children. This is not a toy. Adult supervision is required. * Use only water or solutions recommended by your dentist. * Clean the device and its accessories regularly according to the manufacturer's instructions. * Unplug the device after use and before cleaning. * ****In Case of Problems:** ** Stop using the device immediately and consult a dental health professional.

This document is generated automatically in response to European product safety regulations and in accordance with the GDPR.

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