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****English****

****Safety Warning: ** This oral irrigator is an oral hygiene device. Please read this manual carefully before use.**

*** **Intended Use: ** This device is intended for household use for oral hygiene purposes. Do not use for other purposes. * **At-Risk Individuals: ** Children should be supervised when using this device. Individuals with sensitive gums or dental problems should consult a dentist before use. * **Safety Instructions: ** * Do not immerse the device in water. * Use only water or oral solutions recommended by your dentist. * Do not direct the water jet towards the eyes or ears. * Regularly check the power cord and water tank for any damage. * Unplug the device after use and before cleaning. * Do not modify or repair the device yourself. Contact a qualified professional. * **Cleaning and Maintenance: ** Clean the water tank and nozzle regularly after each use. Use a soft, damp cloth.**

This document is automatically generated in response to European product safety regulations and in compliance with the GDPR.

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