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****ENGLISH****

****Important Safety Warning****

Before using this anti-cellulite device, please read these instructions carefully. Failure to follow these instructions may result in skin irritation, injury, or reduced treatment effectiveness.

* **Contraindications:** Do not use this device if you are pregnant, have blood circulation problems, varicose veins, skin diseases, or wear a pacemaker. Consult a doctor before use if you have any doubts. * **Usage:** Use the device only on clean, dry skin. Do not use on irritated, burned, or injured skin. Do not apply excessive pressure. * **Duration and Frequency:** Do not exceed the recommended usage time in the user manual. Excessive use can cause irritation. * **Maintenance:** Clean the device after each use according to the manufacturer's instructions. Do not immerse the device in water. * **Storage:** Keep the device out of reach of children. * **Possible side effects:** Redness, slight bruising, or skin sensitivity may occur. If these effects persist or worsen, discontinue use and consult a doctor.

This document is automatically generated in response to European product safety regulations and in accordance with the GDPR.

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